



**FRESENIUS
MEDICAL CARE**

**2008[®] SERIES
TECHNICAL TRAINING**

STUDENT GUIDE AND WORKBOOK

The 2008K, 2008K², 2008K@Home™, & 2008T machines constitute the 2008[®] Series Hemodialysis Machines. In this manual the 2008T, 2008K, 2008K² & 2008K@Home™ hemodialysis machines may be referred to as “2008 Series”.

WELCOME TO THE FRESENIUS MEDICAL CARE TRAINING CENTER

i. Typical Daily Class schedule:

- a. 8:30 AM - class starts
- b. 12:00-12:30 - Lunch break
- c. 5:00 PM - class ends
- d. Breaks - two ten-minute breaks in the morning; one in the afternoon

ii. Ground Rules:

- a. Keep your cellphone in silent mode or turn it off
- b. Be punctual to class and from breaks
- c. No exposed legs or open-toe footwear (cut-off jeans and flip-flops)
- d. No web-surfing or personal emails during class

iii. Reference Documents:

- a. 490114 2008K Level 1 Training Manual
- b. 508990-01 Training Presentation: Introduction to Hemodialysis
- c. 508990-02 Training Presentation: Hydraulics- Primary Side
- d. 508990-03 Training Presentation: Hydraulics- Secondary Side
- e. 507297 2008K Preventive Maintenance Procedures
- f. 507781 2008K Semi-Annual Preventive Maintenance Procedures

- iv. Training material handed out in class:
 - a. Hydraulic flow path diagrams
 - b. Calibration resistors for temperature sensor
 - c. Potentiometer adjusting tool
- v. **CLASSROOM COPIES ONLY – PLEASE DO NOT WRITE ON OR TAKE HOME THE FOLLOWING:**
 - a. Level one Training manual
 - b. Calibration procedures
 - c. Preventive Maintenance procedures.
 - d. Description of debug screens
 - e. *biBag*® Technician Manual
- vi. **For a copy of the classroom ONLY documents and class agenda, please refer to Fresenius Medical Care website:**

<https://fmcna.com/product-support-documentation/>

I. Introduction to Dialysis

Learning Objectives:

1. _____
2. _____
3. _____
4. _____

II. Kidney Functions:

1. _____
2. _____
3. _____
4. _____
5. _____

III. Hemodialysis Basics:

The HEMODIALYSIS MACHINE performs four major functions:

1. PREPARES AND MONITORS _____
2. MONITORS THE PATIENT _____
3. CONTROLS _____ FLOW
4. REMOVES _____ (ULTRAFILTRATE) FROM THE PATIENT'S BLOOD
 1. Dialyzer _____
 - i. Blood Compartment _____

 - ii. Dialysate Compartment _____

2. Dialysate Preparation

3. Blood Circuit

4. Controls Dialysate Circuit Flow

The _____ consist of two cylindrical chambers, each with volume 30 (+/- 1)ml
_____ #20 loads one of the balancing chamber with **FRESH** dialysate
_____ #21 loads one of the balancing chamber with **SPENT** dialysate
_____ #22 removes water from patient circuit 1ml per stroke of the pump

5. Diffusion: The passage of _____ from an area of high concentration to an area of low concentration. _____

6. Osmosis: The passage of **WATER** through a semi permeable membrane from an area with a _____ **pressure** gradient to one with a _____ **pressure** gradient.

7. KUF _____

8. 2008 Series Blood Circuit Modules:

Arterial Blood Pump: _____

Level Detector:

Venous Pressure : _____


Air detector: _____

Venous Line Clamp: _____

Optical detector: _____


Heparin Pump: _____

11. Major Machine Sub-Assemblies:




I. MACHINE OPERATION

- a. Machine Function - The 2008 Series is designed to provide hemodialysis treatment by controlling and monitoring both the dialysate and extracorporeal blood circuits
- b. Front View:




FOR TRAINING PURPOSES ONLY – NOT INTENDED FOR USE IN SERVICING MACHINE

c. Rear View:



FEATURE	2008T	2008K²	2008K	2008K@Home
Touchscreen	X		X	X
Keyboard	X			
Touchpad	X			
CDX capability	X			
On-Line Clearance	X	X	X	
Dialysate on/off key		X	X	
Bypass led		X	X	
WI-FI	X			
Ethernet port	X			
USB ports	X			
RS232 port	X	X	X	X

HYDRAULIC BLOCK DIAGRAM



II. HYDRAULIC ASSEMBLY: PRIMARY SIDE

Main function of Primary side: to make _____

- a. Inlet Water Regulator #61:

Manually calibrated to approximately _____ when incoming water valve #41 is closed.

- b. Heat Exchanger #77:

The heated 'spent' fluid, on it way to the drain passes through the exchanger and warms the plate which _____ as it is passing through.

c. Valve #41(27) (incoming water valve):

Opens when the float bob (#5, in hydrochamber C) _____; Closes when the float bob _____.

d. Solenoid valve characteristics:

1. Controlled by the _____
2. Valves are normally _____
3. Energized with ___ volts DC, valve opens
4. Approximately ___ Ohms.

e. Hydro block/Hydrochamber (five chambers, A, B, C, D, E):

f. Heater Element #54:

1300 Watt, _____ element with an internal **resistance** of approximately _____ **ohms**. Voltage is switched to the heater by **a triac** (located in the power supply).

g. NTC #2 (Heater Control Sensor):

Measures temperature and _____ #54 via the triac

h. Float Switch #5:

controls the **water level in the** _____.

i. Deaeration #20:

A strong gear-type pump magnetically coupled to a DC motor controlled by the _____ board.

j. Restrictor Orifice #48:

Restricts _____ to the deaeration pump #20 in Dialysis Program.

k. Loading Pressure Valve #65:

Located at the bottom side of chamber A this valve is manually calibrated to open when the deaeration pump's #20 output pressure reaches **approximately _____ psi.**

l. Concentrate (acid) Pump #16 and Bicarbonate Pump #17

These are piston-diaphragm pumps that must draw and _____ to accurately _____ profiles

m. Acid Filter #71 and Bicarbonate Filter #72

Located inside the _____. They prevent dirt from being drawn into the pump that may damage internal soft parts

n. Mixing Chambers #82

The _____ and _____ concentrates are injected into the loading pressure _____ stream, post hydroblock/hydrochamber and then stirred by the **mixing** chambers. This is the final step of _____ before the Balancing Chambers

- o. #11 Reed Switch Concentrate Wand and # 12 Reed Switch Bicarbonate Wand
Hall effect (reed) switch located at the _____ and _____ rinse port
-

III. HYDRAULIC ASSEMBLY SECONDARY SIDE

Main function of the Secondary side: to deliver _____

- a. Balancing Chamber #68:

The balancing chambers consist of two cylindrical **30 ml** (approximately) chambers. Each chamber contains a non-permeable, _____ that provides complete isolation between _____ **(‘F’, pre-dialyzer) and** _____ **(‘S’ post-dialyzer) fluid**

Valve Cycle 1: Valves # ____, ____, ____, ____ open; Valves # ____, ____, ____, ____ closed.

Valve Cycle 2: Valves # ____, ____, ____, ____ closed ; Valves # ____, ____, ____, ____ open .

Balancing Chamber Valve Dead Time: _____

- b. CFS Pressure Transducer #10:

Strain gauge transducer notifies the Actuator board that it is time to switch the Balancing Chamber valves into the _____ . _____

- c. Conductivity Cell #7 (pre-dialyzer):

The _____ of the fluid flowing through the cell affects the AC resistance path between the probes which in turn varies AC frequency. Frequency is converted to _____ **(mS/cm)** which is displayed digitally to the front panel.

d. NTC #3 (Temperature Monitor Sensor, pre-dialyzer):

Changes in temperature (i.e. changes in voltage) are converted to standard temperature units (**degrees Celsius**) and then displayed digitally to the front panel. NTC#3 values are also _____ readings from Conductivity Cell #7.

e. Heater Runaway Protection:

In the event of a heater control runaway when temperature increases to _____ (way beyond the target value) the _____ (located in the power supply) switches off.

f. Diasafe® PLUS filter Test Valve #28:

is used to _____ the Diasafe® PLUS filter by allowing air into Diasafe® PLUS filter through filter #92.

g. Diasafe Filter #90:

Removes _____ **and** _____ material prior to entering the dialyzer

h. Valve #24 (Dialyzer valve) and Valve #26 (Bypass valve):

Valve# _____ open if _____ **and** _____ readings are within limits

Valve# _____ open if _____ or _____ reading are NOT within limits

i. External Flow Indicator #75 and External Filter #73:

When the machine is _____ (i.e. valve #24 open, #26 closed), dialysate enters the dialyzer which causes the 'bob' to rise and fall.

j. Shunt Door Interlock Switches:

If the shunt door is open ' _____ ' (i.e. valve #24 closes, #26 opens) and also valve #25 closes

k. External Filter #73:

_____ from re-entering the hydraulics via the dialyzer

l. Valve #25:

Solenoid valve that remains open unless the machine is performing

_____ (PHT) or the _____

m. Dialysate Pressure Transducer #9:

This strain gauge transducer senses pressure in the _____.

TMP= _____ - _____

n. Blood Leak Detector #8:

Optically monitors for _____ that may leak through the dialyzer membrane.

o. Post-Dialyzer Temperature Monitor NTC #44:

This monitors temperature post dialyzer. Its only purpose is to _____ readings from Conductivity Cell #13.

p. Post-Dialyzer Conductivity Cell #13:

Used by the machine's _____ function,

q. Level sensor #6:

The “_____” is a hydraulic function that removes air from the spent dialysate before it enters the balancing chamber

r. Air Removal Chamber #69:

Captures and _____ in the fluid prior to entering the Balancing Chambers

s. Valve #43:

Solenoid valve that opens during '_____' and at various times in the Cleaning Programs to disinfect/rinse it flow path.

t. Flow Pump #21:

is a _____ that is magnetically coupled to a DC motor controlled by the Actuator-Test board that precisely controls _____ through the dialyzer

u. CFS Pressure Transducer #10:

v. #78 Flow Pump Relief Valve:

This valve prevents the flow pump from creating _____. It is mechanically **adjusted to _____ PSI or _____ PSI** (depending upon if the machine is Diasafe® PLUS Filter equipped).

w. Valve #30 (drain valve):

The **drain** valve allows spent dialysate to enter _____.

x. Pre-UF Pump Filter #74:

y. Ultrafiltration (UF) Pump #22: The pump is mechanically calibrated to deliver _____/stroke. It is responsible for _____ (weight) from the patient

z. UF Pump Check valve #63 and #64:

The purpose of check valve # ____ is to isolate the UF Pump from the sample port.


The purpose of check valve # ____ is to prevent back flow from the sample port

aa. Valve #39: Opens only in the _____. When open restrictor orifice #48 is bypassed and allows deaeration pump #20 to draw fluid directly out of chamber C. The result is a drastic decrease in deaeration pressure (drops to between 0 and -15 inHg). _____

a.

bb. Valve #29 (recirculation valve):

Allows fluid to the drain to be _____; used in heat disinfect to heat up faster and to conserve energy



a. Start Up Screen

b. Dialysis Screen

c. The Control Panel

d. Treatment Display Section

e. Pressure and Alarm Test

f. Online Clearance

g. Concentrate Select Menu

h. UF Profile Menu

i. Blood Pump Module

j. Level Detector Module

IV. Electronic Circuit Board Review

Warning! Electro Static Discharge PRECAUTIONS REQUIRED

a. Power Supply

- i. Incoming power is stepped down to 24 Volts DC
- ii. Heater Triac
- iii. Power Cord
 - 1. Connects machine to 120V AC power source
 - 2. Provides 120V AC to Heater connections in distribution board
- iv. Power Control Board
 - 1. 12 Volt Standby
 - 2. Heater Relay
 - 3. 24 Volt Relay
 - 4. Mains Fuses (6.3 Amp)

- b. Power Logic Board
 - i. Machine Power On Signal From Power Button
 - ii. Modified versions (now with -12V inverter incorporated)
 - iii. Auto Off option
 - iv. 9 Volt Battery Test
 - v. DC to DC Converters (+5 & +12 Volt DC “regulators”)
 - vi. Heater Control signals
 - vii. Alarm Tones Generated

- c. Actuator & Actuator / Test combination Board
 - i. Drives Pumps, Valves & Motors for Hydraulic control
 - ii. Has Its Own **Analog to Digital Converter (ADC)**
 - iii. **ADC** for communication with Functional board
 - iv. Own software
 - 1. Utilizes Flash Memory
 - 2. Software version not compatible between machine models.
 - v. Secondary Monitor for Hardware Alarm Limits
 - vi. Actuator / Test Combo board **TEST** circuit
 - 1. Used during POST (Power On Self Test)
 - 2. Runs Alarms Test

- d. Functional Board
 - i. Controls all machine Functions (MPU)
 - ii. Watchdog Circuit
 - iii. Calibration EEPROM
 - iv. LCD Drivers
 - v. Software Upgrades are performed using a PAL2 device.
 - vi. All Options are set in Service Mode
 - vii. Runs Pressure Holding Tests
 - 1. Negative Pressure
 - 2. Positive Pressure

e. Test Board

NOTE will not be present if machine equipped with Actuator / Test Combo Board (see above)

- i. Used during POST (Power On Self Test)
- ii. Runs Alarms Test
- iii. EPROM software

f. Sensor Board

- i. Primary Monitor Board
- ii. Has Primary **Analog** to **Digital Converter** for communication with the Functional Board

g. UI / MICS / CDX board -**T machine ONLY!**

- i. UI = User Interface / MICS = Medical Information Computer System / CDX = Clinical Data Exchange
- ii. Transfers screen data from the functional board and up-scales it to the 15 inch LCD display.
- iii. UI portion provides interface between user input and electronic control (Functional board)
- iv. MICS & CDX portion incorporates vendor software installed on the internal flash drive to manage medical information collected or entered during treatment

h. Motherboard

NOTE Motherboards manufactured post September 2007 have additional pin connections (3rd row pins connected @ Act. /Test combo board).

- i. Connections for and communication between circuit boards in card cage.
- ii. Remote connections for Modules.
- iii. Main 24 Volt Power cable connection.
 - 1. Must be plugged in completely!

i. 12 Volts Inverter

NOTE: This board may be incorporated into the Power Logic board (see above)

i. Used by Colin Module

1. NIBP

ii. Used by 2008K machines with OLC

j. Display assembly 2008T

i. LCD graphic display

ii. Touch Screen

iii. Keyboard

iv. Keypad

v. Touch pad

vi. Status Light

k. Display assembly 2008K

i. LCD Display

ii. Display interface board

iii. Front Panel Switch Matrix

iv. Touch Screen

l. Display assembly 2008K²

i. LCD Display

ii. Display interface board

iii. Front Panel Assembly with additional switch matrix buttons

1. + / - Keys

2. Directional arrow keys around 'CONFIRM' button

CLASS EXERCISES AND LABS

Formulas and examples:

1. $TMP = P_{VEN} - P_{DIAL}$
2. $UFR \text{ (ml/Hr)} = [*Pre\text{-}Weight - Dry\ Weight] \div \text{Dialysis Time}$
3. $KUF = \text{ml/Hr}/TMP \text{ (mmHg)}$
4. $KUF = UFR/TMP = UFR \div TMP$
5. $UFR = KUF \times TMP$
6. $TMP = UFR/KUF = UFR \div KUF$
7. $1 \text{ Kg} = 1 \text{ Liter} = 1000 \text{ ml} = 2.2 \text{ lbs.}$

EXAMPLE 1.0

A male patient comes in for dialysis weighing 57 Kg. His dry weight is 53 Kg. He receives no saline and will dialyze for 4 hours on a dialyzer that has a KUF of 5 ml/Hr/mmHg. His P_{VEN} is 160 mmHg. Determine:

- 1) What is his pre-dialysis weight in pounds?
- 2) How much weight in Kg must be removed during the treatment?
- 3) What UFR is required to accomplish this?
- 4) What will the TMP be?
- 5) What will the dialysate pressure be (P_{DIAL})?

ANSWERS:

- 1) $57 \text{ Kg} \times 2.2 \text{ lbs./1 Kg} = \underline{125.4 \text{ lbs.}}$
- 2) $PW - DW = (57 - 53) \text{ Kg} = \underline{4 \text{ Kg}}$
- 3) $UFR = PW - DW \div DT = (57 - 53)\text{Kg} \div 4 \text{ Hr} = 1 \text{ Kg/Hr} = \underline{1000 \text{ ml/Hr}}$
- 4) $TMP = UFR \div KUF = 1000 \div 5 = \underline{200 \text{ mmHg}}$
- 5) Rearranging Equation 2.0, $P_{DIAL} = P_{VEN} - TMP$ **Equation 2.6**

$$P_{DIAL} = P_{VEN} - TMP = 160 - 200 = \underline{-40 \text{ mmHg}}$$

EXAMPLE 1.1

If the dialyzer, from example 1.0, was changed to one that has a KUF of 50 calculate the: 1) TMP; 2) Dialysate Pressure

ANSWERS:

1) $TMP = UFR \div KUF = 1000 \div 50 = \underline{20 \text{ mmHg}}$

2) $P_{DIAL} = P_{VEN} - TMP = 160 - 20 = \underline{+140 \text{ mmHg}}$

Exercise:

EXAMPLE ONLY – NOT TO BE USED IN SERVICING MACHINE

Calibration & Identification Check List

IDENTIFICATIONS:

- Electronics card cage ID complete *initials* _____
- Main Power supply ID complete..... *initials* _____
- Inlet Water Regulator #61 *initials* _____
- #77 Heat Exchanger *initials* _____
- Solenoid Valve #41 / *initials* _____
- Solenoid Valve #29..... *initials* _____
- Chamber A *initials* _____
- Heater & NTC#2 *initials* _____
- The Float in Chamber C..... *initials* _____
- Bypass Valve for Deaeration *initials* _____
- Concentrate Pumps *initials* _____
- Filters for Acid & Bicarb intake *initials* _____
- Acid & Bicarb Rinse Ports & Reed Switches *initials* _____
- Mixing Chamber or Chambers *initials* _____
- Balancing Chamber Assembly..... *initials* _____
- Conductivity Cell#7 & NTC#3 *initials* _____
- Diasafe Filter & Diasafe Test Valve *initials* _____
- Solenoid Valve #24..... *initials* _____
- Solenoid Valve #26..... *initials* _____
- External Flow Indicator *initials* _____
- Shunt Door and Switches..... *initials* _____
- External Line Filter *initials* _____
- Valve 25..... *initials* _____
- Dialysate Pressure Transducer..... *initials* _____
- Blood Leak Detector..... *initials* _____
- Post Conductivity#13 and NTC#44..... *initials* _____
- Air Separation Chamber..... *initials* _____
- Flow Pump *initials* _____
- Chamber Full Switch (CFS) *initials* _____
- Flow Relief Valve *initials* _____
- Drain Valve *initials* _____
- Pre – UF Filter..... *initials* _____
- UF Pump..... *initials* _____
- Check Valves *initials* _____

EXAMPLE ONLY – NOT TO BE USED IN SERVICING MACHINE

Calibration & Identification Check List (Student's Copy) Cont....

Hydraulic calibrations

1. Inlet Water Regulator - Inlet water pressure 18 - 20psi (**dialysis mode**) **Calibration complete**
Enter reading from pressure gauge _____ *initials* _____
2. Deaeration & Loading pressure
 - a. Deaeration ***at sea level*** -24inHg
Enter reading from pressure gauge _____ *initials* _____
 - b. Loading Pressure 23-25psi
18-20psi without Diasafe
Enter reading from pressure gauge _____ *initials* _____
3. Flow Pump Relief Pressure 35 – 36 psi
29-30psi without Diasafe
Enter reading from pressure gauge _____ *initials* _____
4. Balance Chamber volume
Enter volume measured _____ *initials* _____
5. Acid pump volume
Enter volumes measured 1) _____ 2) _____ *initials* _____
6. Bicarbonate pump volume
Enter volumes measured 1) _____ 2) _____ *initials* _____
7. UF Pump volume 1cc/stroke
Enter volume measured _____ *initials* _____

Sensor Calibrations

1. Arterial Pressure Calibration *initials* _____
2. Venous Pressure Calibration *initials* _____
3. Dialysate Pressure Calibration *initials* _____
4. Temperature Sensor Calibration *initials* _____
5. Post Temperature Sensor Calibration (if applicable) *initials* _____
6. Temperature Control Calibration *initials* _____
7. Blood Leak Calibration *initials* _____
8. Conductivity Cell(s) Calibration *initials* _____

Monitor Calibrations

1. Set Clock (service mode or dialysis under B.P. screen) *initials* _____
2. Voltage Detection *voltage reading* _____ *initials* _____
3. Arterial Pump Rate *initials* _____

CALIBRATION DONE IN DIALYSIS MODE!

- Level Detector Calibration (channel 1 and channel 2) (**dialysis mode**) *initials* _____