Optiflux®

F180NR

Hollow Fiber Dialyzers

NOT FOR REUSE

GENERAL INFORMATION

Indications: Optiflux F180NR dialyzers are designed for single use acute and chronic hemodialysis.

USA only: Federal law restricts this device to sale by or on order of a physician.

Contraindications: Specific contraindications for the dialyzer are unknown. Generally, the contraindications for hemodialysis are applicable. The dialyzer should only be used as directed by a physician.

Precautions: In the event of a blood leak during dialysis, the health care provider should respond according to the facility's established protocol.

Air entering the extracorporeal circuit during dialysis can result in serious injury or death. Check the integrity of all blood lines prior to the initiation of dialysis and periodically during the treatment. The venous return line or drip chamber should be continuously monitored with an air detector.

Warning: Due to the high water flux capability of high permeability membranes with a KUF ≥ 8, it is necessary to use such dialyzers only in conjunction with dialysis machines that are equipped with precise ultrafiltration control, such as the Fresenius A2008 series. We recommend that dialyzers with a KUF ≥ 6 only be used with such UF control machines. In any case, the safety instructions of the manufacturer of the hemodialysis machine must be followed.

The user is cautioned to regularly monitor the patient's chemistry values using quantitative measurements and analysis to ensure that the expected therapy is delivered. The clinical parameters monitored should, at least, include: urea, hematocrit, and serum albumin.

Dialyzers may leak, resulting in patient blood loss or contamination with dialysate.

Dialysate: The dialysate must meet AAMI standards for dialysis (RD5). The use of bicarbonate-containing dialysate is recommended for short dialysis.

Side effects: In rare cases hypersensitivity reactions may occur during hemodialysis treatment. A history of allergies is an indication for careful monitoring of hypersensitivity reactions. Dialyzers of this type should not be used again on any patient exhibiting a hypersensitivity reaction. With severe reactions, dialysis must be discontinued and aggressive, first line therapy for anaphylactoid reactions must be initiated. The decision to return the blood must be made by a physician.

Heparinization: It is recommended to systemically heparinize the patient by allowing the prescribed loading dose of heparin to circulate for 3 to 5 minutes before beginning extracorporeal circulation. In addition, the extracorporeal circuit may be heparinized, taking care that the total amount of heparin to be delivered meets the prescription. During dialysis, the dosage of heparin and the mode of administration are the responsibility of the attending physician. The coagulation time should be checked regularly.

Sterile/Non-pyrogenic: The dialyzer blood pathway is sterile and non-pyrogenic if the blood port caps are in place and undamaged. Sterilized with ethylene oxide gas. Do not use if the dialyzer is damaged in any way. Use aseptic technique for all blood side connections. Structural integrity of the hemodialyzer is warranted for the first use only when prepared as directed.

Recommended storage: Between 5 and 30°C (41 and 86°F).

Dialyzer reuse: Optiflux F180NR dialyzers are **NOT** designed for reuse.

PREPARATION FOR DIALYSIS - DRY PACK

- Place the dry dialyzer in a vertical position, arterial end downward.
- Install the arterial and venous bloodlines on the hemodialysis machine. Note: Refer to the manufacturer's instructions for the dialysate delivery machine.
- Remove any dialyzer blood port caps and aseptically connect the arterial and venous blood lines to the dialyzer.
- Aseptically spike a 1 liter bag of 0.9% sterile normal saline with a clamped IV administration set. Attach the IV administration set to the patient end of the arterial bloodline.
- Open the clamp on the IV set. Prime the arterial bloodline, dialyzer, and venous bloodline using a blood pump speed of approximately 150 mL/min. Discard the first 500 mL of solution. The drip chambers should be maintained about 3/4 full.
- Stop the blood pump. Clamp the arterial and venous bloodlines. **Turn the dialyzer so that the venous end is downward.** Aseptically connect the patient ends of the arterial and venous lines together in preparation for recirculation. Open the clamps on the bloodlines.
- Verify that the dialysate is within the prescribed conductivity limits with a calibrated external conductivity meter. To identify situations where the acetate or acid and bicarbonate concentrates are not properly matched, use pH paper or a meter to verify that the approximate pH is in the physiologic range.
- Attach the dialysate lines to the dialyzer. Fill the dialysate compartment. In order to maximize the efficiency of the dialyzer, the dialysate flow must be countercurrent to the blood flow.
- Rotate the dialyzer so that the arterial end is downward. Recirculate the blood side at a flow rate of 300- 400 mL/min and a dialysate flow of 500 mL/min for a minimum of ten to fifteen minutes and until all the air has been purged from the system before connecting to the patient. Continue recirculation and dialysate flow until patient connection.
- Ultrafilter or flush an additional 500 cc of 0.9% sterile normal saline so that the extracorporeal circuit has been flushed with a minimum of 1 L of saline to minimize sterilization residues.

If the dialysate delivery system was chemically disinfected or sterilized prior to patient use, be sure to test for the absence of germicide residuals with a test intended for this application, according to the test manufacturer's instructions.

INITIATION OF DIALYSIS

- Turn the blood pump off. Clamp the saline line and the arterial and venous bloodlines.
- Aseptically attach the arterial bloodline to the patient's arterial access. Open the arterial bloodline clamp.
- Place the venous bloodline in a drain container, making sure not to contaminate the end of the bloodline. Open the clamp on the venous bloodline.
Warning: This step must be carefully monitored to prevent any possibility of blood loss.
- Turn the blood pump speed up to 100 - 150 mL/min. Fill the extracorporeal circuit with the patient's blood.
- Do not infuse the recirculated saline prime into the patient. If saline is required for volume enhancement, discard the recirculated saline and fill the bloodlines with fresh saline.
- Turn off the blood pump and clamp the venous bloodline.
- Aseptically attach the patient end of the venous bloodline to the patient's venous access.
- Open the clamp to the venous access.
- Unclamp the venous bloodline and set the blood flow to the prescribed rate. Rotate the dialyzer so that the venous end is downward.
- Set the prescribed ultrafiltration rate. A minimum ultrafiltration rate of 300 mL/hr during dialysis treatment is recommended with Optiflux dialyzers.

DURING DIALYSIS TREATMENT:

- If a blood leak should occur during the treatment, the decision to attempt to allow the leak to clot off by reducing the blood flow and ultrafiltration rate to minimum values is a clinical decision. The decision whether or not to return the blood to the patient must be made by a medical professional.
- Air entering the extracorporeal circuit during dialysis is very serious and should be avoided. A routine check of all connections prior to initiation of dialysis and periodically throughout the treatment is recommended. Constant monitoring of the venous return line with an air detector is essential. Should air get into the venous line during treatment, the dialysis must be discontinued without returning any of the patient's blood that is mixed with air.

TERMINATION OF DIALYSIS

- When the dialysis treatment is completed, turn the blood flow rate to zero and UF rate to recommended minimum.
- Clamp arterial bloodline and arterial access, and aseptically disconnect from the patient's arterial access.
- Using the blood pump, rinse the patient's blood back using sterile 0.9% saline solution at a slow rate. **Do not allow air to enter the extracorporeal circuit.**
- Once the blood has been returned, turn the blood pump flow rate to zero.
- Clamp the venous bloodline.
- Clamp the patient's venous access and aseptically disconnect the venous bloodline from the patient's access.
- Discard the extracorporeal circuit in an appropriate biohazard waste receptacle.** References: 29CFR 1910.145, 1910.1030 (Code of Federal Regulations) and appropriate state or local codes.

TECHNICAL DATA: These data represent typical *in vitro* performance. Actual *in vivo* performance will differ.

	F180NR	
Ultrafiltration coefficient (<i>in vitro</i> , bovine blood, Hct 32%)	55 413	mL/h/mmHg mL/h/kPa
Clearance (Q _b = 300 mL/min, Q _d = 500 mL/min)		
Urea*	270	mL/min
Creatinine	245	mL/min
Phosphate	238	mL/min
Vitamin B ₁₂	165	mL/min
Priming volume blood	105	mL
Flow resistance blood (Q _b = 200 mL/min)	47 6.3	mmHg kPa
Flow resistance - dialysis fluid (Q _d = 500 mL/min)	16 2.1	mmHg kPa
Maximum TMP	600 80	mmHg kPa
Maximum blood flow	600	mL/min
Maximum dialysate flow	1000	mL/min
Surface area	1.8	m ²

*Sodium used as a marker for Urea.

Note: Clearance tests performed using aqueous solutions of sodium, creatinine, phosphate, and Vitamin B₁₂.

Membrane material:	Advanced Fresenius Polysulfone®
Inner diameter (nominal):	200 µm
Housing:	Polycarbonate
Potting compound:	Polyurethane
O-ring:	Silicone
Blood connectors:	DIN 13090 Part 3
Dialysis fluid connectors:	DIN 58352 Part 2

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