# Optiflux®

## F18NR F16NR

### Hollow Fiber Dialyzers

### NOT FOR REUSE

#### GENERAL INFORMATION

**Indications:** Optiflux F18NR and F16NR dialyzers are designed for single use acute and chronic hemodialysis.

**USA only:** Federal law restricts this device to sale by or on order of a physician.

**Contraindications:** Specific contraindications for the dialyzer are unknown. Generally, the contraindications for hemodialysis are applicable. The dialyzer should only be used as directed by a physician.

**Precautions:** In the event of a blood leak during dialysis, the health care provider should respond according to the facility's established protocol.

Air entering the extracorporeal circuit during dialysis can result in serious injury or death. Check the integrity of all blood lines prior to the initiation of dialysis and periodically during the treatment. The venous return line or drip chamber should be continuously monitored with an air detector.

**Warning: Due to the high water flux capability of high permeability membranes with a KUF > 8, it is necessary to use such dialyzers only in conjunction with dialysis machines that are equipped with precise ultrafiltration control, such as the Fresenius A2008 series. We recommend that dialyzers with a KUF > 6 only be used with such UF control machines. In any case, the safety instructions of the manufacturer of the hemodialysis machine must be followed.**

The user is cautioned to regularly monitor the patient's chemistry values using quantitative measurements and analysis to ensure that the expected therapy is delivered. The clinical parameters monitored should, at least, include: urea, hematocrit, and serum albumin.

Dialyzers may leak, resulting in patient blood loss or contamination with dialysate.

**Dialysate:** The dialysate must meet AAMI standards for dialysis (RD5). The use of bicarbonate-containing dialysate is recommended for short dialysis.

**Side effects:** Adverse reactions such as hypertension, hypotension, headache, and nausea may be associated with hypovolemia or hypervolemia and can usually be avoided with careful management of patient fluid, electrolyte balance, blood flow rate, and ultrafiltration rate.

In rare cases, thrombocytopenia or hypersensitivity reactions including anaphylactic or anaphylactoid reactions to the dialyzer, or other elements in the extracorporeal circuit, may occur during hemodialysis. Hypersensitivity reactions may cause mild to severe signs and symptoms, including: itching, flushing, hives, swelling, fever, leukopenia, hypotension, hypertension, shortness of breath with wheezing, arrhythmias, and/or respiratory arrest.

Patients with a history of hypersensitivity reactions or patients who have a history of being highly sensitive and allergic to a variety of substances should be carefully monitored during treatment. With severe hypersensitivity reactions, dialysis must be discontinued and aggressive first line therapy for hypersensitivity reactions must be initiated. The decision to return the patient's blood in the event of a hypersensitivity reaction is determined by the physician.

**Heparinization:** It is recommended to systemically heparinize the patient by allowing the prescribed loading dose of heparin to circulate for 3 to 5 minutes before beginning extracorporeal circulation. During dialysis, the dosage of heparin and the mode of administration are the responsibility of the attending physician. The coagulation time should be checked regularly.

**Sterile/Non-pyrogenic:** The dialyzer blood pathway is sterile and non-pyrogenic if the blood port caps are in place and undamaged. Electron beam sterilized. Do not use if the dialyzer is damaged in any way. Use aseptic technique for all blood side connections. Structural integrity of the hemodialyzer is warranted for the first use only when prepared as directed.

**Recommended storage:** Between 5 and 30°C (41 and 86°F).

**Dialyzer reuse:** Optiflux F18NR and F16NR dialyzers are **NOT** designed for reuse.

#### PREPARATION FOR DIALYSIS - DRY PACK

- Place the dry dialyzer in a vertical position, arterial end downward.
- Install the arterial and venous bloodlines on the hemodialysis machine. Note: Refer to the manufacturer's instructions for the dialysate delivery machine.
- Remove any dialyzer blood port caps and aseptically connect the arterial and venous blood lines to the dialyzer.
- Aseptically spike a 1 liter bag of 0.9% sterile normal saline with a clamped IV administration set. Attach the IV administration set to the patient end of the arterial bloodline.
- Open the clamp on the IV set. Prime the arterial bloodline, dialyzer, and venous bloodline using a blood pump speed of approximately 150 mL/min. Discard the first 500 mL of solution. The drip chambers should be maintained about 3/4 full.
- Stop the blood pump. Clamp the arterial and venous bloodlines. **Turn the dialyzer so that the venous end is downward.** Aseptically connect the patient ends of the arterial and venous lines together in preparation for recirculation. Open the clamps on the bloodlines.
- Verify that the dialysate is within the prescribed conductivity limits with a calibrated external conductivity meter. To identify situations where the acetate or acid and bicarbonate concentrates are not properly matched, use pH paper or a meter to verify that the approximate pH is in the physiologic range.
- Attach the dialysate lines to the dialyzer. Fill the dialysate compartment. In order to maximize the efficiency of the dialyzer, the dialysate flow must be countercurrent to the blood flow.
- Rotate the dialyzer so that the arterial end is downward. Recirculate the blood side at a flow rate of 300- 400 mL/min and a dialysate flow of 500 mL/min for a minimum of ten to fifteen minutes and until all the air has been purged from the system before connecting to the patient. Continue recirculation and dialysate flow until patient connection.

- Ultrafilter or flush an additional 500 cc of 0.9% sterile normal saline so that the extracorporeal circuit has been flushed with a minimum of 1 L of saline to minimize sterilization residues.
- If the dialysate delivery system was chemically disinfected or sterilized prior to patient use, be sure to test for the absence of germicide residuals with a test intended for this application, according to the test manufacturer's instructions.
- Do not infuse the recirculated saline prime into the patient. Discard the recirculated saline and fill the entire extracorporeal circuit with fresh saline prior to connecting to the patient. The volume of fresh saline used to fill the extracorporeal circuit should be equal to the volume of the dialyzer and blood tubing set in use.

#### INITIATION OF DIALYSIS

- To initiate dialysis; stop the blood pump, clamp the dialysis priming set and the arterial and venous bloodlines.
- Aseptically attach the patient ends of the bloodlines to the patient's arterial and venous access. Open the arterial and venous bloodline clamps and the clamps on the patient access.
- Increase the blood pump speed slowly to the prescribed blood flow rate. Be sure to monitor the arterial and venous blood pressures carefully during this process to note any possible flow restrictions or inappropriate pressure readings.
- Once the prescribed blood flow rate has been achieved, set the prescribed ultrafiltration rate.

#### DURING DIALYSIS TREATMENT

- If a blood leak should occur during the treatment, the decision to attempt to allow the leak to clot off by reducing the blood flow and ultrafiltration rate to minimum values is a clinical decision. The decision whether or not to return the blood to the patient must be made by a medical professional.
- Air entering the extracorporeal circuit during dialysis is very serious and should be avoided. A routine check of all connections prior to initiation of dialysis and periodically throughout the treatment is recommended. Constant monitoring of the venous return line with an air detector is essential. Should air get into the venous line during treatment, the dialysis must be discontinued without returning any of the patient's blood that is mixed with air.

#### TERMINATION OF DIALYSIS

- When the dialysis treatment is completed, turn the blood flow rate to zero and UF rate to recommended minimum.
- Using a plastic clamp, clamp the arterial blood tubing just above the saline "T" connection (between the saline "T" connection and the inlet of the blood pump segment).
- Open the clamps on the saline line and saline "T" connection and allow the 0.9% saline solution to rinse the blood between the saline "T" and the patient arterial access connection back to the patient.
- Move the clamp from above the saline "T" connection to below the saline "T" connection and clamp the patient's arterial access and close the main line clamp on the arterial bloodline.
- Using the blood pump, slowly return the blood in the rest of the extracorporeal circuit to the patient.
- Do not allow air to enter the extracorporeal circuit.**
- Once the blood has been returned, turn the blood pump flow rate to zero.
- Clamp the main line clamp on the venous bloodline.
- Clamp the patient's venous access.
- Aseptically disconnect the arterial and venous bloodlines from the patient's access.
- Discard the extracorporeal circuit in compliance with appropriate federal, state and local requirements, which may require the use of medical waste and/or biohazard receptacles. References: 29CFR, 1910.145, 1910.1030 (Code of Federal Regulations) and appropriate state or local codes

**TECHNICAL DATA:** These data represent typical *in vitro* performance. Actual *in vivo* performance will differ.

	F18NR	F16NR	
Ultrafiltration coefficient ( <i>in vitro</i> , bovine blood, Hct 32%)	12 90	10 75	mL/h/mmHg mL/h/kPa
Clearance (Qb = 300 mL/min, Qd = 500 mL/min)			
Urea*	265	262	mL/min
Creatinine	236	228	mL/min
Phosphate	188	175	mL/min
Vitamin B <sub>12</sub>	113	100	mL/min
Priming volume blood	103	84	mL
Flow resistance blood (Qb = 200 mL/min)	53 7.1	56 7.5	mmHg kPa
Flow resistance - dialysis fluid (Qd = 500 mL/min)	11 1.5	14 1.9	mmHg kPa
Maximum TMP	600 80	600 80	mmHg kPa
Maximum blood flow	600	500	mL/min
Maximum dialysate flow	1000	1000	mL/min
Surface area	1.8	1.5	m <sup>2</sup>

\*Sodium used as a marker for Urea.

Note: Clearance tests performed using aqueous solutions of sodium, creatinine, phosphate, and

Vitamin B<sub>12</sub>

Membrane material: Advanced Fresenius Polysulfone®

Inner diameter (nominal): 200 µm

Housing: Polycarbonate

Potting compound: Polyurethane

O-ring: Silicone

Blood connectors: DIN 13090 Part 3

Dialysis fluid connectors: DIN 58352 Part 2

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