Instructions for Use
OPTIFLUX®
F160NR F180NR F200NR F250NR
Capillary High Flux Dialyzer
Single Use Only - Electron Beam Sterilized


Figure 1: Dialyzer Orientation During Blood Compartment Priming

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
Figure 2: Dialyzer Orientation for Hemodialysis Treatment with Arterial-end up.

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
Figure 3: Dialyzer Orientation for Hemodialysis Treatment with Venous-end up.

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GENERAL INFORMATION

Indications For Use: Optiflux Dialyzers are intended for patients with acute or chronic renal failure when conservative therapy is judged to be inadequate.

Contraindications: Optiflux dialyzers are not intended for use in patients that weigh less than 30kg. **USA Only:** Federal (US) law restricts this device to sale by or on order of a physician.

Caution: The operator should strictly adhere to the manufacturer's recommended procedures, warnings, and cautions as listed in these instructions for use.

Precautions: Dialyzers may leak resulting in patient blood loss or contamination with dialysate. In the event of a blood leak during dialysis, the health care provider should respond according to the facility's established protocol.

Air entering the extracorporeal circuit during dialysis can result in serious injury or death. Check the security of all extracorporeal connections prior to the initiation of dialysis and periodically throughout the treatment. The venous drip chamber should be continuously monitored with a level / air detector.

Warning: Due to the high water flux capability of high permeability membranes with an ultrafiltration coefficient > 8, it is necessary to use such dialyzers only in conjunction with dialysis machines that are equipped with precise ultrafiltration control, such as the Fresenius 2008 series. In any case, the safety instructions for the hemodialysis machine must be followed.

Dialyzer performance may change with the duration of observation.

The user is cautioned to regularly monitor the patient's chemistry values using quantitative measurements and analysis to ensure that the prescribed therapy is delivered. The clinical parameters monitored should at least include urea, hemoglobin and serum albumin.

Dialysate: The dialysate must meet current standards for dialysis AAMI RDS2 and ISO 11663.

Side effects: Adverse reactions such as hypertension, hypotension, headache, and nausea may be associated with hypovolemia or hypervolemia and can usually be avoided with careful management of patient fluid, electrolyte balance, blood flow rate, and ultrafiltration rate.

In rare cases, thrombocytopenia or hypersensitivity reactions including anaphylactic or anaphylactoid reactions to the dialyzer, or other elements in the extracorporeal circuit, may occur during hemodialysis. Hypersensitivity reactions may cause mild to severe signs and symptoms, including: itching, flushing, hives, swelling, fever, leukopenia, hypotension, hypertension, shortness of breath with wheezing, arrhythmias, and/or respiratory arrest.

Patients with a history of hypersensitivity reactions or patients who have a history of being highly sensitive and allergic to a variety of substances should be carefully monitored during treatment. With severe hypersensitivity reactions, dialysis must be discontinued and aggressive first line therapy for hypersensitivity reactions must be initiated. The decision to return the patient's blood in the event of a hypersensitivity reaction is determined by the physician.

Heparinization: We recommend systemic heparinization of the patient. Systemic heparinization is the administration of the prescribed loading dose of heparin into the vascular access of the patient; then waiting a period of 3 to 5 minutes prior to initiating the treatment. During dialysis, the dose of heparin and method of administration is the decision of the physician.

Sterile/Non-pyrogenic: The blood pathway is sterile and non-pyrogenic if blood port caps are in place and undamaged. Do not use if the dialyzer is damaged in any way. Use aseptic technique for all blood side connections. Structural integrity of the hemodialyzer is warranted for the first use only when prepared as directed. **Sterilized using the electron beam method of sterilization.**

Recommended storage: Between 5 and 30 degrees C (41 - 86 degrees F).

Dialyzer reuse: Optiflux dialyzers are not designed for or intended for reuse.

PREPARATION FOR DIALYSIS

• If the dialysate delivery system was chemically disinfected or sterilized prior to patient use, be sure to test for the absence of germicide residuals with a test intended for this application, according to the test manufacturer's instructions.

• Place the dialyzer in the dialyzer holder in the vertical position.

• Install the arterial and venous bloodlines on the hemodialysis machine.

Note: Refer to dialysate delivery machine manufacturer's instructions for use for setting up bloodlines.

• Remove the blood port cap from the arterial end of the dialyzer and aseptically connect the arterial (red) dialyzer end of the bloodline to the lower dialyzer blood port.

• Remove the blood port cap from the venous end of the dialyzer and aseptically connect the venous (blue) dialyzer end of the bloodline to the upper dialyzer blood port.

• Ensure all connections are secure. Refer to Figure 1: Position of Dialyzer During Blood Compartment Priming.

• Aseptically spike a 1 liter bag of 0.9% sterile saline solution with a clamped dialysis priming set.

• If not already attached, attach the dialysis priming set to the saline "T" connection located just before the blood pump segment on the arterial bloodline. Check to be sure the connection is secure.

• Open the clamp on the dialysis priming set and allow saline to gravity prime the portion of the arterial bloodline from the saline "T" to the patient end.

• Clamp the main line tubing on the arterial bloodline between the patient end and the saline "T" connection.

• Start the blood pump and set a pump speed of 150mL/min. Prime the remainder of the arterial bloodline, dialyzer, and venous bloodline with saline. While the extracorporeal circuit is filling with saline, intermittently pinch and release the bloodline between the blood pump and the dialyzer to help purge air from the dialyzer. Gently tap the dialyzer to facilitate air removal from the dialyzer.

• Fill the dialyzer and bloodlines with 300mL sterile 0.9% saline solution. The drip chambers in the bloodlines should be set to and maintained at ¾ full.

• Stop the blood pump. Clamp the venous bloodline. Aseptically connect the patient ends of the arterial and venous bloodlines together in preparation for recirculation of the extracorporeal circuit. Unclamp main line clamps on arterial and venous bloodlines.

• Perform Pressure and Alarm tests on Fresenius 2008 machine.

• Verify that the dialysate is within the prescribed conductivity limits with a calibrated independent conductivity meter. To identify situations where the acetate or acid and bicarbonate concentrates are not properly matched, use a calibrated pH meter to verify that the pH of the dialysate is within the appropriate physiologic range.

• Attach the dialysate lines to the dialysate ports of the dialyzer. In order to maximize the efficiency of the dialyzer assure that the dialysate lines are attached so that the blood and dialysate flows are counter-current.

• Rotate the dialyzer so the venous end is lower than the arterial end to fill the dialysate compartment with dialysate.

• Assure that the dialysate compartment is completely filled with dialysate.

• Return the dialyzer to the arterial end down position.

• Recirculate the extracorporeal circuit at a blood flow rate of 300 to 400mL/min and a dialysate flow 500mL/min until all air has been purged from the dialyzer and bloodlines.

• During recirculation, to assist in removing air from the dialyzer, intermittently pinch and release the blood tubing between the blood pump and the dialyzer. It is recommended to gently tap the upper portion of the dialyzer to remove air even if the header area of the dialyzer looks free of air.

• Do not infuse the recirculated saline prime into the patient. Discard the recirculated saline and fill the entire extracorporeal circuit with fresh saline prior to connecting to the patient. The volume of fresh saline used to fill the extracorporeal circuit should be equal to the volume of the dialyzer and blood tubing set in use.

• If a blood leak should occur during the treatment, the operator should follow established facility procedures for a dialyzer blood leak.

• Aseptically attach the patient ends of the bloodlines to the arterial and venous access of the patient. Open venous mainline and vascular access clamps.

• Turn on the blood pump and slowly increase the blood pump speed to the prescribed blood flow rate. Be sure to monitor the hemodialysis machine arterial and venous blood pressures carefully during this process to note any possible flow restrictions or inappropriate pressure readings.

• Once the prescribed blood flow rate has been achieved, set the prescribed ultrafiltration rate.

• Performance testing has demonstrated that the dialyzer can be run in either the arterial end up or arterial end down position during the hemodialysis treatment. Refer to Figure 2 or 3: Hemodialysis Treatment for proper dialyzer and bloodline positioning.

DURING THE DIALYSIS TREATMENT

• If a blood leak should occur during the treatment, the operator should follow established facility procedures for a dialyzer blood leak.

• Air entering the extracorporeal circuit during dialysis is a very serious event and should be avoided. A routine check of all connections prior to initiation of dialysis and periodically throughout the dialysis treatment is recommended. Constant monitoring of the venous drip chamber with a level / air detector is required. Should air enter the venous line during the treatment, the dialysis treatment must be discontinued without returning any of the blood mixed with air.

COMPLETION OF DIALYSIS

• When the dialysis treatment is completed, turn the blood pump off and set the UF rate to the recommended minimum. Verify sufficient 0.9% sterile saline solution for return of the blood in the extracorporeal circuit back to the patient.

• Open the priming line clamp and place a hemostat on the arterial bloodline between the saline "T" and the blood pump. Return the blood in the mainline between the saline "T" and the patient end back to the patient.

• Close the arterial mainline clamp when return is completed.

• Remove the hemostat and start the blood pump at 100-150 mL/min pump speed. Intermittently pinch and release the blood tubing between the blood pump and the dialyzer to help to efficiently return the blood in the extracorporeal circuit back to the patient. Do not let air enter the extracorporeal circuit.

• Once the blood has been returned to the patient, turn the blood pump off. Clamp the venous mainline and the patient's arterial and venous access. Aseptically disconnect the arterial and venous bloodlines from the patient's access.

• Discard the extracorporeal circuit in an appropriate biohazard waste receptacle. References: 29CFR, 1910.145, 1910.1030 (Code of Federal Regulations) and appropriate state and local codes.

Technical Data: These data represent typical *in vitro* performance. Actual *in vivo* performance may differ.

| | Optiflux F160NR | Optiflux F180NR | Optiflux F200NR | Optiflux F250NR |
|--|-----------------|-----------------|-----------------|-----------------|
| Catalogue Number | 0500316E | 0500318E | 0500320E | 0500325E |
| Membrane Surface Area (m²) | 1.5 | 1.7 | 1.9 | 2.5 |
| Ultrafiltration Coefficient¹ (mL/hr/mmHg), Qb=300mL/min | 61 | 76 | 74 | 111 |
| Priming Volume Blood (mL) | 87 | 102 | 113 | 142 |
| KoA (mL/min) Qb=300mL/min, Qd=500 mL/min | 1167 | 1321 | 1415 | 1714 |
| Flow Resistance (mmHg) | | | | |
| Blood Qb=300mL/min | 89 | 80 | 66 | 57 |
| Dialysate Qd=500mL/min | 17 | 15 | 15 | 19 |
| Maximum TMP (mmHg) | 600 | 600 | 600 | 600 |
| Minimum Blood Flow (mL/min) | 100 | 100 | 100 | 100 |
| Maximum Blood Flow (mL/min) | 600 | 600 | 600 | 600 |
| Minimum Dialysate Flow (mL/min) | 200 | 200 | 200 | 200 |
| Maximum Dialysate Flow (mL/min) | 800 | 800 | 800 | 800 |
| Maximum Ultrafiltration Rate (mL/hr) | 4000 | 4000 | 4000 | 4000 |

| Clearance² (mL/min) | | | | |
|---|-----|-----|-------|-----|
| Qb=300mL/min, Qd=500mL/min, | | | | |
| Urea ³ | 271 | 277 | 280 | 287 |
| Creatinine | 242 | 253 | 259 | 273 |
| Phosphate | 241 | 253 | 254 | 274 |
| Vitamin B ₁₂ | 153 | 169 | 177 | 205 |
| Lysozyme ⁴ | 74 | 87 | 83 | 113 |
| Qb=400mL/min, Qd=500mL/min, | | | | |
| Urea ³ | 317 | 328 | 333 | 349 |
| Creatinine | 274 | 290 | 298 | 321 |
| Phosphate | 279 | 297 | 298 | 328 |
| Vitamin B ₁₂ | 162 | 181 | 189 | 226 |
| Qb=400mL/min, Qd=600mL/min, | | | | |
| Urea ³ | 332 | 344 | 348 | 366 |
| Creatinine | 288 | 304 | 310 | 338 |
| Phosphate | 290 | 307 | 309 | 341 |
| Vitamin B ₁₂ | 167 | 189 | 195 | 235 |
| Qb=400mL/min, Qd=800mL/min, | | | | |
| Urea ³ | 349 | 360 | 365 | 381 |
| Creatinine | 302 | 320 | 328 | 357 |
| Phosphate | 304 | 321 | 324 | 356 |
| Vitamin B ₁₂ | 173 | 195 | 203 | 248 |
| Sieving Coefficient⁵ Qb=300mL/min, ultrafiltration rate=29mL/min | | | | |
| Albumin | | | <0.01 | |

Dialyzer Technical Information

| | |
|---|---------------------------------|
| Membrane Material: | Advanced Fresenius Polysulfone® |
| Fiber Inner Diameter: | 185 microns |
| Membrane Wall Thickness -- F160NR, F180NR, F200NR: | 35 microns |
| Membrane Wall Thickness -- F250NR: | 30 microns |
| Housing Material: | Polycarbonate |
| Potting Compound: | Polyurethane |
| O-ring: | Silicone |
| Recommended Blood Connections: | ISO 8637:2010 |
| Recommended Dialysis Fluid Connections: | ISO 8637:2010 |
| Sterilization Method: | Electron Beam |

Footnotes:

¹ Bovine blood, Hct 32%, protein 60g/L, 37°C

² Aqueous solution of respective test substance, 37°C; ultrafiltration rate=0mL/min

³ Sodium used as a marker for urea

⁴ Lysozyme used as a marker for middle molecule solutes; molecular weight 14,300 Daltons

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