

Fresenius Medical Care North America Service Bulletin

Equipment: 2008T Hemodialysis Machines
 Bulletin: 19-FRT-001 Rev B
 Subject: 2008T Machines With New Style Hydraulics

1.0 PURPOSE

To inform all Field Service and Technical Support personnel as well as Fresenius Medical Care-trained customers of the use of a new style hydraulics in 2008[®]T hemodialysis machines.


2.0 DESCRIPTION

2008T machines are being manufactured with new style hydraulics. Either the new or current style hydraulics may be present in a 2008T machine as we transition through inventory. The table below shows the starting serial number sequence for each machine part number which are equipped with this new style hydraulics. The new style hydraulics can also be identified by their part numbers (F40005844 and F40005846).

Machine Part Number	Beginning Serial Number
190713	8T0S233715
190766	8T0S226306
190858	8T0S226703
190895	8T0S226464

The functionality of these hydraulics are identical to the previous with the following differences:

- The acid and bicarbonate inlet lines have been extended to 8 feet in length to ease in the use of a central concentrate system. If a central concentrate system is not used, these lines can be shortened to the previous length.
- These hydraulics include extended life components inside the acid, bicarbonate and Ultrafiltration pumps. These components qualify the machine for a 2 year preventive maintenance (PM) schedule.
- The valves have been upgraded from an “EKS” valve to a “Ratio” valve. They can be identified in the figures below. These valves can be used interchangeably in the hydraulics.



If there are any questions regarding this bulletin, contact Fresenius Medical Care Technical Support at 800-227-2572.

© 2019, Fresenius Medical Care North America, All rights reserved.