Hemoflow™

F3, F4

Hollow Fiber Dialyzers

NOT FOR REUSE

GENERAL INFORMATION

Indications: F3, F4 dialyzers are designed for single use in acute and chronic hemodialysis.

USA only: Federal law restricts this device to sale by or on order of a physician.

CAUTION: The operator should strictly adhere to the manufacturer's recommended procedures, warnings, and cautions as listed in these instructions for use.

Contraindications: Specific contraindications are unknown. Generally the contraindications for hemodialysis are applicable. The dialyzer should only be used as directed by a physician.

Precautions: Dialyzers may leak resulting in patient blood loss or contamination with dialysate. In the event of a blood leak during dialysis, the health care provider should respond according to the facility's established protocol.

Air entering the extracorporeal circuit during dialysis can result in serious injury or death. Check the security of all extracorporeal connections prior to the initiation of dialysis and periodically throughout the treatment. The venous drip chamber should be continuously monitored with a level detector.

Warning:

- The user is cautioned to regularly monitor the patient's chemistry values using quantitative measurements and analysis to ensure that the prescribed therapy is delivered. The clinical parameters monitored should at least include urea, hemoglobin, and serum albumin.
- Use only an ultrafiltration-controlled hemodialysis machine such as the Fresenius 2008 series.

Dialysate: The dialysate must meet AAMI standards for dialysis (RD 52).

Side effects: In rare cases, hypersensitivity reactions to the dialyzer or other elements in the extracorporeal circuit may occur during hemodialysis. If a hypersensitivity reaction occurs, the source of the hypersensitivity should be identified and that component of the extracorporeal circuit should be excluded from future use in hemodialysis treatments for that patient. With severe reactions, dialysis must be discontinued and aggressive first line therapy for hypersensitivity reactions must be initiated. The decision to return the patient's blood in the event of a hypersensitivity reaction is the decision of the physician.

Anticoagulation: Unless medically contraindicated, it is recommended to anticoagulate the extracorporeal circuit during hemodialysis. The choice of anticoagulant, the dose, and method of administration is the decision of the physician.

Sterile/Non-pyrogenic: The dialyzers are sterilized using ethylene oxide gas. The dialyzer blood pathway is sterile and non-pyrogenic if the blood port caps are in place and undamaged. Do not use if the dialyzer is damaged in any way. Use aseptic technique for all bloodside connections. Structural integrity of the dialyzer is warranted for the first use only when prepared as directed.

Recommended storage: Between 5 and 30° C. (41 – 86° F).

Dialyzer reuse: Hemoflow F3, F4 dialyzers are not designed or intended for reuse.

PREPARATION FOR DIALYSIS – DRY PACK

- Place the dialyzer in the dialyzer holder in the vertical position, arterial end downward. Install the appropriate arterial and venous bloodlines on the hemodialysis machine.
Note: Refer to dialysate delivery machine manufacturer's instructions for use in setting up bloodlines.
- Remove blood port caps from the dialyzer and aseptically connect the arterial and venous dialyzer ends of the bloodlines to the dialyzer. Check to be sure connections are secure.
- Aseptically spike a 1 liter bag of 0.9% sterile saline solution with a clamped dialysis priming set.
- If not already attached, attach the dialysis priming set to the saline "T" connection, located just before the blood pump segment on the arterial bloodline. Check to be sure the connection is secure.
- Open the clamp on the dialysis priming set and allow saline to gravity prime the portion of the arterial bloodline from the saline "T" to the patient end.
- Clamp the main line tubing on the arterial bloodline between the patient end and the saline "T" connection.
- Start the blood pump and set a pump speed of 150 mL/min. Prime the rest of the arterial bloodline, dialyzer and venous bloodline with saline. While the extracorporeal circuit is filling with saline, intermittently pinch and release the bloodline between the blood pump and the dialyzer to help to facilitate removal of air from the dialyzer and extracorporeal circuit.
- Discard the first 500 mL of saline solution. The bloodline drip chamber levels should be set to and maintained at ¾ full.
- Stop the blood pump. Clamp the arterial and venous bloodlines. Aseptically connect the patient ends of the arterial and venous bloodlines together in preparation for recirculation of the extracorporeal circuit. Unclamp main line clamps on arterial and venous bloodlines.
- Perform Pressure Holding Test to check the integrity of the dialysis equipment's UF control system if using the Fresenius 2008 or other volumetric ultrafiltration control machine.
- Verify that the dialysate is within the prescribed conductivity limits with a calibrated independent conductivity monitor. To identify situations where the acetate or acid and bicarbonate concentrates are not properly matched, use a calibrated pH meter to verify that the pH of the dialysate is within the appropriate physiologic range.
- Rotate the dialyzer so the venous end is down. Attach the dialysate lines to the dialyzer. Fill the dialysate compartment with the dialyzer in the venous end down position. In order to maximize the efficiency of the dialyzer, the dialysate flow must be counter-current to the blood flow.
- When the dialysate compartment is filled, turn the dialyzer back to the arterial end down position and place back in dialyzer holder.
- Recirculate the extracorporeal circuit at a blood flow rate of 300 mL/min and a dialysate flow of 500 mL/min. Continue to recirculate the circuit for 10 to 15 minutes and until all visible air has been purged from the dialyzer and bloodlines before connecting to the patient. Continue recirculation and dialysate flow until patient connection.
- To assist in removing air from the dialyzer during recirculation, intermittently pinch and release the blood tubing between the blood pump and the dialyzer. Even if the header area of the dialyzer looks free of air, it is recommended that gentle tapping of the venous end of the dialyzer, with the arterial end down, should be done to remove air from the dialyzer.
- Ultrafilter and/or flush an additional 500 mL of 0.9% sterile normal saline so that the extracorporeal circuit has been flushed with a minimum of 1 liter of saline to minimize sterilization residues.
- If the dialysate delivery system was chemically disinfected or sterilized prior to patient use, be sure to test the dialysis machine for the absence of germicide residuals with a test intended for this application, according to the test manufacturer's instructions.
- Do not infuse the recirculated saline prime into the patient. Discard the recirculated saline and fill the bloodlines with fresh saline prior to connecting to the patient. The volume of fresh saline used to fill the extracorporeal circuit should be equal to the volume of the dialyzer and bloodline set in use.

INITIATION OF DIALYSIS

- To initiate dialysis: stop the blood pump, clamp the dialysis priming set, and the arterial and venous bloodlines.
- Aseptically attach the patient ends of the bloodlines to the patient's arterial and venous access. Open the arterial and venous bloodline clamps and the clamps on the patient access.
- Increase the blood pump speed slowly to the prescribed blood flow rate. Be sure to monitor the arterial and venous blood pressures carefully during this process to note any possible flow restrictions or inappropriate pressure readings.
- Once the prescribed blood flow rate has been achieved, set the prescribed ultrafiltration rate.

DURING THE DIALYSIS TREATMENT

- If a blood leak should occur during the treatment, the operator should follow the facility's established procedure for a dialyzer blood leak.
- Air entering the extracorporeal circuit during dialysis is a very serious event and should be avoided. A routine check of all connections prior to initiation of dialysis and periodically throughout the dialysis treatment is recommended.
- Constant monitoring of the venous drip chamber with a level detector is required. Should air get into the venous line during the treatment, the dialysis must be discontinued without returning any of the blood mixed with air.
- Continuously monitor fluid removal using an independent device such as a crib or bed scale.

TERMINATION OF DIALYSIS

- When the dialysis treatment is completed, turn the blood pump off and set the UF rate to the recommended minimum. Check to see that there is enough 0.9% saline in the saline bag for rinsing the blood in the extracorporeal circuit back to the patient.
- Using a hemostat, clamp the arterial bloodline between the saline "T" and the blood pump. Rinse the blood in the tubing between the saline "T" and the patient end back to the patient.
- Clamp the arterial bloodline between the patient connection and the saline "T". Remove the clamp on the bloodline between the saline "T" and the blood pump.
- Start the blood pump and set at a 150 to 200 mL/min pump speed. Intermittently pinch and release the blood tubing between the blood pump and the dialyzer to help to efficiently rinse the blood in the extracorporeal circuit back to the patient. Do not let air enter the extracorporeal circuit during rinseback.
- Once the blood has been returned to the patient, turn the blood pump off. Clamp the arterial and venous bloodlines and the patient's arterial and venous access. Aseptically disconnect the arterial and venous bloodlines from the patient's access.
- Discard the extracorporeal circuit in an appropriate biohazard waste receptacle. References: 29CFR, 1910.145, 1910.1030 (Code of Federal Regulations) and appropriate state and local codes.

Technical data: These data represent typical *in vitro* performance. Actual *in vivo* performance may differ.

	F3	F4	
Ultrafiltration coefficient (<i>in vitro</i> , bovine blood, Hct=32%)	1.7	2.8	mL/hr/ mmHg
Clearance Q _b =200mL/min, Q _d =500mL/min UF=0mL/min			
Urea*	125	155	mL/min
Creatinine	95	128	mL/min
Phosphate	47	76	mL/min
Vitamin B ₁₂	20	32	mL/min
Priming volume blood	24	42	mL
Flow resistance blood (Q _b = 200 mL/min)	178	97	mmHg
Flow resistance dialysate (Q _d = 500 mL/min)	18	21	mmHg
Maximum TMP	600	600	mmHg
Maximum blood flow	200	300	mL/min
Maximum dialysate flow	800	800	mL/min
Surface area	0.4	0.7	m ²

*Sodium used as marker for urea.

Note: Clearance tests performed using aqueous solution of sodium, creatinine, phosphate, and Vitamin B₁₂

Membrane material	Fresenius Polysulfone®
Fiber inner diameter (nominal)	200 microns
Membrane wall thickness	40 microns
Housing	Polycarbonate
Potting compound	Polyurethane
O-ring	Silicone
Blood connections	DIN 13090 Part 3
Dialysis fluid connections	DIN 58352 Part 2
Sterilization Method	Ethylene Oxide

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