Instructions for Use
OPTIFLUX®
F180NR
High Flux Dialyzer
Single Use Only — Ethylene Oxide Sterilized


Figure 1: Dialyzer Orientation During Blood Compartment Priming.

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
Figure 2: Dialyzer Orientation for Hemodialysis Treatment with Arterial-end up.

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
Figure 3: Dialyzer Orientation for Hemodialysis Treatment with Venous-end up.

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Instructions for Use
OPTIFLUX®
F180NR
Capillary High Flux Dialyzer
Single Use Only — Ethylene Oxide Sterilized

GENERAL INFORMATION

Indications For Use: Optiflux Dialyzers are intended for patients with acute kidney injury or chronic kidney disease when conservative therapy is judged to be inadequate.

Contraindications: Specific contraindications are unknown. Generally the contraindications for hemodialysis are applicable. The dialyzer should only be used as directed by a physician. Patients with known hypersensitivity to any of the materials used (see technical data) must not be treated with this product.

CAUTION: Federal law restricts this device to sale by or on the order of a physician. Selection of the appropriate dialyzer should take into consideration underlying patient health conditions. Caution for use of the dialyzer should be considered in certain patients, including patients with low weight or blood volume, who are elderly, a history of hypotension, who are pregnant, with severe heart failure and, on life support.

Provide increased caution if using the device in patients recently initiating hemodialysis.

To protect the dialyzer integrity, it is necessary to avoid direct exposure to sunlight during storage.

Instructions for Use: The operator should strictly adhere to the manufacturer's recommended procedures, warnings, and cautions as listed in these instructions for use.

Precution: Dialyzers may leak resulting in patient blood loss or contamination with dialysate. In the event of a blood leak during dialysis, the health care provider should respond according to the facility's established protocol.

Warning:

Air entering the extracorporeal circuit during the dialysis treatment phase can result in air embolism, serious injury or death. Check the security of all extracorporeal connections prior to the initiation of dialysis and at regular intervals throughout the treatment. The venous drip chamber should be continuously monitored with a level/air detector. In the event that air enters the extracorporeal circuit and cannot be removed safely, discontinue the treatment. In such a circumstance, do not return the blood to the patient.

The user is cautioned to regularly monitor the patient's chemistry values using quantitative measurements and analysis to ensure that the prescribed therapy is delivered. The clinical parameters monitored should at least include urea, hemoglobin and serum albumin.

Use only with dialysis machines that are equipped with precise ultrafiltration control due to the high permeability membranes (ultrafiltration coefficient >8) of these dialyzers. Failure to do so could lead to hypovolemia. Only use the device in dialysis machines that have intact blood leak detectors and air detectors.

Do not use a dialyzer past expiration date. Doing so may result in compromised dialyzer performance, foreign body embolism, or compromised sterility.

Do not use a dialyzer with visible damage, missing or unattached blood port caps to ensure patient safety. Doing so may result in compromised sterility or blood loss.

Do not infuse recirculated saline into the patient. Doing so may result in toxicity to the patient.

Initiate treatment only after dialysate compartment is completely filled with dialysate. Air pockets in the dialysate fluid compartment may lead to insufficient clearance.

If a patient is receiving oral or intravenous medications during a treatment, refer to the medication label for instructions on dosage.

Dialysate: The dialysate must meet current water quality standards for dialysis (e.g. ISO 11663) to mitigate the risk of infection due to back-filtration.

Side effects: Adverse reactions such as hypertension, hypotension, headache, and nausea may be associated with hypovolemia or hypervolemia and can usually be avoided with careful management of patient fluid, electrolyte balance, blood flow rate, and ultrafiltration rate. Hemolysis may also occur as a result of dialysis.

In rare cases, thrombocytopenia or hypersensitivity reactions including anaphylactic or anaphylactoid reactions to the dialyzer, or other elements in the extracorporeal circuit, may occur during hemodialysis. Hypersensitivity reactions may cause mild to severe signs and symptoms, including: itching, flushing, hives, swelling, fever, leukopenia, hypotension, hypertension, shortness of breath with wheezing, arrhythmias, and/or respiratory arrest.

Patients with a history of hypersensitivity reactions or patients who have a history of being highly sensitive and allergic to a variety of substances should be carefully monitored during treatment. With severe hypersensitivity reactions, dialysis must be discontinued and aggressive first line therapy for hypersensitivity reactions must be initiated. Blood from the extracorporeal system must not be returned to the patient in cases of hypersensitivity reactions.

Heparinization: We recommend systemic heparinization of the patient. Systemic heparinization is the administration of the prescribed loading dose of heparin into the vascular access of the patient; then waiting a period of 3 to 5 minutes prior to initiating the treatment. During dialysis, the dose of heparin and method of administration is the decision of the physician.

Sterile/Non-pyrogenic: The blood pathway is sterile and non-pyrogenic if blood port caps are in place and undamaged. Do not use if the dialyzer is damaged in any way. Use aseptic technique for all blood side connections. Structural integrity of the hemodialyzer is warranted for the first use only when prepared as directed. **Sterilized using the ethylene oxide method of sterilization.**

Recommended storage: Between 5 and 27 degrees C (41 – 80 degrees F).

Dialyzer reuse: Optiflux dialyzers are single use only. Re-use of the device may create a potential risk to the patient, including contamination and/or impairment of the device function.

Note: Follow the hemodialysis machine specific operating instructions for treatment, when applicable; otherwise, refer to the instructions below.

PREPARATION FOR DIALYSIS

- Verify dialyzer selected matches with treatment order.
- Remove dialyzer from packaging.
- Attach dialysate port caps to dialyzer.
Note: Do not remove blood port caps until ready to attach arterial and venous bloodlines to the dialyzer.
- Place the dialyzer in the dialyzer holder in the vertical position.
- Install the arterial and venous bloodlines on the hemodialysis machine.
Note: Refer to dialysate delivery machine manufacturer's instructions for use for setting up bloodlines.
- Remove bottom blood port cap from the dialyzer and discard.
- Aseptically connect the arterial bloodline to the dialyzer.
- Remove remaining blood port cap from the dialyzer and discard.
- Aseptically connect the venous bloodline to the dialyzer.
- Refer to Figure 1: Dialyzer Orientation During Blood Compartment Priming.
- Ensure all connections are tight.
- Hang 1 liter bag of 0.9% sterile saline solution and spike with the saline line.

Filling Blood Compartment

- If not already attached, attach the dialysis priming set to the saline "T" connection located just before the blood pump segment on the arterial bloodline. Check to be sure the connection is secure.
- Allow saline to gravity prime the arterial bloodline.
- Close the clamp on the arterial bloodline when complete.
- Start the blood pump and set a pump speed of 150mL/min to prime the remainder of the arterial bloodline, dialyzer, and venous bloodline with saline.
Note: While the extracorporeal circuit is priming with saline, gently tap the dialyzer and intermittently pinch and release the bloodline to help purge air from the dialyzer.
- Prime and discard 500 mL of 0.9% saline solution through the dialyzer. The bloodline drip chamber levels should be set to and maintained at ¾ full.
- Stop the blood pump after discarding 500 mL of saline.
- Aseptically connect the patient ends of the arterial and venous bloodlines together in preparation for recirculation of the extracorporeal circuit.
- Open arterial and venous bloodline clamps.

Filling Dialysate Fluid Compartment

- Remove the dialysate port caps from the dialyzer and save.
- Attach the dialysate lines to the dialysate ports of the dialyzer. In order to maximize the efficiency of the dialyzer, assure that the dialysate lines are attached so that the blood and dialysate flows are counter-current.

- Rotate the dialyzer to arterial end up to fill the dialysate compartment with dialysate.

Recirculation

- Ensure entire bloodline is ready for recirculation including connections, clamp position, and drip chamber fill volume.
- Recirculate the extracorporeal circuit at a blood flow rate of at least 150 mL/min and a dialysate flow of at least 300 mL/min until all air has been purged from the dialyzer and bloodline.
- During recirculation, intermittently pinch and release the bloodline tubing to assist in removing air from the dialyzer. It is recommended to gently tap the venous end of the dialyzer to remove air even if the header appears to be free of air.
- When ready to initiate treatment, allow saline to gravity flush the arterial bloodline.
- Use the blood pump to flush the venous bloodline with 500 mL of 0.9% sterile saline solution so that a minimum of 1 liter of saline has been primed and flushed through the dialyzer to minimize residues.
Warning: Do not infuse the recirculated saline into the patient. Doing so may result in toxicity to patient.

INITIATION OF DIALYSIS

Warning: Initiate treatment only after dialysate compartment is completely filled with dialysate. Air pockets in the dialysate fluid compartment may lead to insufficient clearance.

- Refer to Figure 2 or Figure 3 for dialyzer orientation during treatment.
- Note: Refer to bloodline manufacturer's instructions for use for appropriate dialyzer orientation during treatment.
- Stop the Blood Pump and close the saline line clamps.
- Aseptically attach the patient ends of the bloodlines to the arterial and venous access of the patient. Open the arterial and venous bloodline and vascular access clamps.
- Turn on the blood pump and slowly increase the blood pump speed to the prescribed blood flow rate. Be sure to monitor the hemodialysis machine arterial and venous pressures carefully during this process to note any possible flow restrictions or inappropriate pressure readings.
- Once the prescribed blood flow rate has been achieved, set the prescribed ultrafiltration rate.

DURING THE DIALYSIS TREATMENT

- If an external blood leak (e.g. loose connection, cracked header/housing) should occur during the treatment, the operator should follow the established facility procedures.
- If an internal blood leak is detected during the treatment, the operator should follow the established facility procedures.

COMPLETION OF DIALYSIS

- When treatment has ended and the blood has been returned to the patient, close the bloodline clamps and aseptically disconnect the arterial and venous bloodlines from the patient's access.
- Remove the dialysate lines from the dialyzer and apply the dialysate port caps.
Note: If dialysate port caps are not available, rotate the dialyzer horizontally, making sure to have the dialysate ports facing up to prevent leaking before removing the dialysate lines from the dialyzer.
- Discard the extracorporeal circuit in an appropriate biohazard waste receptacle.** References: 29CFR, 1910.145, 1910.1030 (Code of Federal Regulations) and appropriate state and local codes.

Technical Data:

These data represent typical *in vitro* performance. Actual *in vivo* performance may differ. Dialyzer performance may change with the duration of observation (treatment).

Catalogue Number	Optiflux F180NR
Membrane Surface Area (m ²)	1.7
Ultrafiltration Coefficient ¹ (mL/hr/mmHg), Q _b =300mL/min	72
Priming Volume Blood (mL)	97
K _o A (mL/min)	1336
Q _b =300mL/min, Q _d =500 mL/min	
Flow Resistance (Pressure Drop) (mmHg)	
Blood Q _b =300mL/min	91
Dialysate Q _d =500mL/min	15
Maximum TMP (mmHg)	600
Minimum Blood Flow (mL/min)	100
Maximum Blood Flow (mL/min)	600
Minimum Dialysate Flow (mL/min)	200
Maximum Dialysate Flow (mL/min)	800
Maximum Ultrafiltration Rate (mL/hr)	4000
Clearance ² (mL/min)	
Q _b =300mL/min, Q _d =500mL/min	
Urea ³	278
Creatinine	247
Phosphate	245
Vitamin B ₁₂	162
Lysozyme ⁴	75
Q _b =400mL/min, Q _d =500mL/min	
Urea ³	321
Creatinine	283
Phosphate	286
Vitamin B ₁₂	174
Q _b =400mL/min, Q _d =600mL/min	
Urea ³	337
Creatinine	296
Phosphate	297
Vitamin B ₁₂	181
Q _b =400mL/min, Q _d =800mL/min	
Urea ³	354
Creatinine	313
Phosphate	309
Vitamin B ₁₂	188
Sieving Coefficient ⁵	
Q _b =300mL/min, ultrafiltration rate=29mL/min	
Albumin	<0.01

Footnotes:

- Bovine blood, Hct 32%, protein 60g/L, 37°C
- Aqueous solution of respective test substance, 37°C; ultrafiltration rate=0mL/min
- Sodium used as a marker for urea
- Lysozyme used as a marker for middle molecule solutes; molecular weight 14,300 Daltons
- Bovine plasma, protein 60g/L, 37°C

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