



# Fresenius Medical Care

## Fresenius Medical Care North America Service Bulletin

Equipment: 2008T Hemodialysis Machines  
 Bulletin: 12-FRT-006 Rev A  
 Subject: CDX & Software Update Test Kit Upgrade

Originated By (Dept.):	Reviewed By	Approved By	Effective Date
J. Fox (TSD)	Kenneth Cogswell	CN#043904	9/12/2012

Distribution: Dir TS (B. La Bier), TS Mgr (K. Cogswell), TS Trng Mgr (E. Cordova), Dir QS (M. Owen), QE (M. Mohr), Eng (C. Woods), QC Mgr (D. McKenzie)  
 Prod Mgr (M. Singer, G. Remy), FMS TS Mgr (J. Barr), Canada (T. Greene), Mexico (J. Naranjo)

### 1.0 PURPOSE

To inform all Field Service and Technical Support personnel of how to upgrade the 2008T CDX & Software Update Test Kit with the latest Linux Operating (OS) system and instructions.

### 2.0 SCOPE

The 2008T CDX & Software Update Test Kit (P/N 190755) can be upgraded without the expense of a new router by using the 2008T CDX Test Kit Upgrade (P/N 190915):

#### 2008T CDX & Software Update Test Kit (P/N 190755\*) includes:

- 190756 2008T USB CDX Test Script 2.0 Flash Drive
- 365075-03 Wireless Router
- 190731 CDX Compact Flash Card with OS
- 508419 2008T CDX & Software Update Test Instructions

#### 2008T CDX Test Kit Upgrade (P/N 190915\*) includes

- 190731 CDX Compact Flash Card with Latest Linux OS
- 508419 2008T CDX & Software Update Test Instructions

**\*Note:** The 2008T CDX & Software Update Test Kit (P/N 190755) is automatically updated to the latest CDX Compact Flash Card with OS. When ordering new test kit, a 2008T CDX Test Kit Upgrade (P/N 190915) is not needed.

© 2012, Fresenius Medical Care North America, All rights reserved.

If there are any questions regarding this bulletin, contact Fresenius Technical Support at 800-227-2572.