



FRESENIUS MEDICAL CARE

Fresenius Medical Care North America Service Bulletin

Equipment: 2008T Hemodialysis Machines
Bulletin: 12-FRT-006 Rev C
Subject: CDX & Software Update Test Kit Upgrade (Update)

1.0 PURPOSE

To inform all Field Service and Technical Support personnel of how to upgrade the 2008T CDX & Software Update Test Kit with the latest Linux Operating (OS) system and instructions.

2.0 SCOPE

The 2008T CDX & Software Update Test Kit (P/N 190755) can be upgraded without the expense of a new router by using the 2008T CDX Test Kit Upgrade (P/N 190915):

2008T CDX & Software Update Test Kit (P/N 190755*) includes:

- 190756 2008T USB CDX Test Script 3.0 Flash Drive
- 365075-03 Wireless Router
- 190731 CDX Compact Flash Card with OS
- 191297 CDX eSATA Solid State Drive with Linux OS
- 508419 2008T CDX & Software Update Test Instructions

2008T CDX Test Kit Upgrade (P/N 190915*) includes

- 190731 CDX Compact Flash Card with Latest Linux OS
- 191297 CDX eSATA Solid State Drive with Latest Linux OS
- 508419 2008T CDX & Software Update Test Instructions

***Note:** The 2008T CDX & Software Update Test Kit (P/N 190755) is automatically updated to the latest CDX Compact Flash Card with OS. When ordering new test kit, a 2008T CDX Test Kit Upgrade (P/N 190915) is not needed.

If there are any questions regarding this bulletin, contact Fresenius Medical Care Technical Support at 800-227-2572.

© 2012 - 2023, Fresenius Medical Care. All Rights Reserved.