



## Fresenius Medical Care North America Service Bulletin

Equipment: Dry Acid Dissolution Units  
Bulletin: 13-DRY-001 Rev A  
Subject: Dry Acid Dissolution Unit Part Numbers

### 1.0 PURPOSE

To inform all Field Service and Technical Support personnel of new part numbers unique to the 99 gallon and 132 gallon Dry Acid Dissolution Unit.

### 2.0 DESCRIPTION

The 99 gallon Dry Acid Dissolution Unit uses many of the same parts as the 99 gallon GranuFlo Dissolution Unit; however, there are some parts that are unique. The table below shows the part number and description of these parts.

PART NUMBER	DESCRIPTION
160136	Assembly, Front Panel Dry Acid, 99 Gal.
160151	Dry Acid 99 System Cable Harness
363029-09	Breaker Power Switch UL
G15-026	Fuse, 250 V 1A, SLO-BLO

The 132 gallon Dry Acid Dissolution Unit uses many of the same parts as the 132 gallon GranuFlo Dissolution Unit; however, there are some parts that are unique. The table below shows the part number and description of these parts.

PART NUMBER	DESCRIPTION
160142	Control Console Assy, Dry Acid 132
160149	Fuse Holder and Fuse Assy
160150	Assembly, Dry Acid 132 Power Cord
363029-09	Breaker Power Switch UL
G15-026	Fuse, 250 V 1A, SLO-BLO

If there are any questions regarding this bulletin, contact Fresenius Medical Care Technical Support at 800-227-2572

© 2013, Fresenius Medical Care North America, All rights reserved.