

Product Information Notice

GranuFlo® Dry Acid Concentrate

March 2020

Dear Valued Customer,

In some GranuFlo® Dry Acid Concentrate cases, it is possible for a small area of orangish discoloration to occur on the surface of the product inside the blend bags. The discoloration is inherent to the manufacturing process, and has been demonstrated to eventually disappear. The contents of the blend bags meet the requirements of the US Pharmacopeia. The safety and efficacy of the GranuFlo product are not affected by this orange discoloration.



Image 1



Image 2



Image 3

NOTE:

- Any type of discoloration noticed in the acid blend bag, apart from that illustrated in **Images 1-3**, should be reported immediately.
- Each case of GranuFlo contains three bags: two acid blend bags and one bag of dextrose



**FRESENIUS
MEDICAL CARE**

RENAL TECHNOLOGIES