



Fresenius Medical Care North America Service Bulletin

Equipment: All 2008 Series Hemodialysis Machines
Bulletin: 16-FRS-003 Rev A
Subject: Extended Lines Option and Dialysate Composition Verification

1.0 PURPOSE

To inform all Field Service and Technical Support personnel as well as Fresenius Medical Care-trained customers of the availability of the following documentation for the 2008® Series hemodialysis machines:

- A procedure on how to install extended acid and bicarbonate concentrate supply lines
- The availability of a dialysate composition verification procedure to be used to verify conductivity after conductivity related repair work is performed.

2.0 SCOPE

The procedure to extend acid and bicarbonate concentrate supply lines is now available for use with all 2008 Series hemodialysis machines.

The procedure to verify dialysate composition after conductivity related repair work is now available for use with all 2008 Series hemodialysis machines.

Please consult the following documents found on the FMCNA.com website:

- 2008 Series Extended Acid & Bicarbonate Tubing Installation Instructions (P/N 490286)
- Dialysate Composition Verification Procedure (P/N 490287)

If there are any questions regarding this bulletin, contact Fresenius Medical Care Technical Support at 800-227-2572.

© 2016, Fresenius Medical Care North America, All rights reserved.